マウス HGF EIA

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Mouse HGF EIA

Mouse HGF determination kit by EIA

- Thoroughly read this instruction manual before use of this kit.
- This kit is for research use only.

I. General Description

Hepatocyte growth factor (HGF), which is a potent growth factor as a trigger for liver regeneration, act as accelerator for growth of various cells not only hepatocyte but also epithelial cells, endothelial cells, and some of mesenchymal tissue. Several studies have reported that HGF had correlation with tumor progression. It also established that HGF enhanced the motility of various types of cells, induced morphogenesis and tubule formation, and controlled the immune response, and many researchers have studied using mice or rats.

This kit is developed for determination of mouse hepatocyte growth factor (mouse HGF) for research purposes. Mouse HGF can be determined quantitatively using the mouse HGF standard solution provided with this kit.

This kit also cross reacts with human HGF but its cross reactivity is about 4%.

II. Kit components


1. Solid phase microplate (8 wells/strip x 12) 1 plate
(Coated with anti-HGF monoclonal antibody)
2. Mouse HGF standard solution 5 vials
(1 vial each of 0, 1, 3, 10 and 25 ng/mL standard concentration: 0.5 mL/vial)
3. Sample diluent L, low salt concentration 15 mL x 1 vial
4. Sample diluent H, high salt concentration 15 mL x 1 vial
5. Anti-HGF rabbit antibody 10 mL x 1 vial
(Anti-HGF rabbit polyclonal antibody)
6. Enzyme labeled antibody 10 mL x 1 vial
(Peroxidase labeled anti-rabbit goat immunoglobulin)
7. Enzyme substrate 10 mL x 1 vial
(3, 3', 5, 5'-tetramethylbenzidine)
8. Reaction stopper 10 mL x 1 vial
9. 20x concentrated washing solution 25 mL x 3 vials
10. Plate seal 4 sheets

III. Application

Determination of mouse hepatocyte growth factor (mouse HGF)

IV. Assay principle

The detection system of this kit is of a sandwich configuration based on enzyme immuno assay consisting of three steps of antigen-antibody reaction between a monoclonal antibody raised against genetically engineered human HGF and most cross-reactive with mouse HGF (anti-HGF mouse monoclonal antibody), mouse HGF in samples, anti-HGF rabbit antibody and anti-rabbit goat immunoglobulin, followed by color development by enzyme reaction. The 1st reaction takes place between the anti-HGF mouse monoclonal antibody immobilized on a solid phase microplate and mouse HGF in samples; the 2nd reaction between mouse HGF in samples bound to the solid plate antibody and anti-HGF rabbit polyclonal antibody, and the third reaction between the enzyme labeled antibody (peroxidase labeled anti-rabbit goat immunoglobulin) and the anti-HGF rabbit polyclonal antibody. After the third reaction, color proportional to the mouse HGF concentration is developed by enzyme reaction and its absorbance is measured for the determination of mouse HGF on a calibration curve prepared using the mouse HGF standard solution.

**V. Operation****1. Preparation of reagent**

- 1) Washing solution
Dilute 20x concentrated washing solution 20 times with purified water. Keep this solution at 2 ~ 10°C after preparation.

2. Materials required but not provided

- 1) One each of micro pipettes, 50 µL and 100 µL
- 2) A measuring cylinder, 1 L
- 3) An aspirator and a polyethylene washing bottle, or a microplate washer
- 4) A dark box (A light tight cupboard or a drawer will do.)
- 5) A microplate reader capable of reading wavelength at 450 nm (reference wavelength at 620 nm or longer)

3. Assay procedure

For the determination of samples^{*1} in high salt concentration buffer, for example, extraction from organs, follow the instruction for high salt concentration samples (**Method A**), and for samples in low salt concentration buffer, for example, plasma^{*2}, follow the instruction for low salt concentration samples (**Method B**). Refer to **Appendix 1** for extraction HGF in the tissue sample, or and **Appendix 2** for purification of it.

Samples which may show absorbance beyond the calibration curve should be diluted with the sample diluent before assay.


*1: Keep samples in plastic vials such as polyethylene and polypropylene. Do not use glass vials.

*2: Platelet is rich with mouse HGF and can be easily damaged during isolation of plasma. To prevent mixture with platelet derived HGF, isolate plasma with absolute care to minimize damage of platelet.

Note:

- Solid phase microplate consists of 12 detachable 8-well strips. Repack in the aluminum pouch unused strips together with desiccant and seal it for storage at 2 ~ 10°C.
- While in use, Mouse HGF standard solution should be maintained at the ice cold water temperature and all reagents should be stored at 2 ~ 10°C after use. Samples should be stored frozen (-20°C or below).
- Return the kit to 15 ~ 30°C before use.
- As illustrated in **Fig.1**, using 2 wells for each standard concentration, measure absorbance of Mouse HGF standard solutions in each assay (2 wells x 5 different standard concentration).

Fig.1 Method A & Method B

**1) Adding samples****(1) High salt concentration samples (Method A)**

If sample dilution is required, dilute with Sample diluent H before assay.

Add 50 µL each of Sample diluent H in the wells for the standard solution (10 wells) and add 50 µL of each Mouse HGF standard solution (concentration 0, 1, 3, 10, and 25 ng/mL), 2 wells for each concentration.

Add 50 µL each of Sample diluent L in sample wells and add 50 µL each of samples, 2 wells for each sample.

(2) Low salt concentration samples (Method B)

If sample dilution is required, dilute with Sample diluent L before assay.

Add 50 µL each of Sample diluent L in all standard solution and sample wells. Add 50 µL of each Mouse HGF standard solutions (concentration 0, 1, 3, 10, and 25 ng/mL) and samples in respective wells, 2 wells for each standard solution and sample.

2) 1st reaction

Cover the microplate with the plate seal provided and leave the microplate to stand at 15 ~ 30°C for 20 hrs (16 ~ 24 hrs).

3) Washing

Remove the plate seal and suck out the well content with an aspirator.

Using a polyethylene washing bottle, fill wells with the washing solution prepared in **1) Washing solution of "1. Preparation of reagent"**. Hold the microplate upside down and vigorously shake out the washing solution. After repeating this washing 5 times, gently tap the plate surface on paper towel to remove the washing solution.

Note: While washing the microplate wells, do not dry the inner wall of wells.

As soon as washing is over, immediately follow the next steps.

4) Addition of Anti-HGF rabbit antibody

Add 100 µL each of Anti-HGF rabbit antibody in all wells.

5) 2nd reaction

Cover the microplate with the plate seal provided. Leave the microplate to stand at 15 ~ 30°C for 2 hrs.

6) Washing


Repeat the procedure **3)** above to wash the microplate wells.

- 7) Addition of Enzyme labeled antibody
Add 100 μ L each of Enzyme labeled antibody in all wells.
- 8) 3rd reaction
Cover the microplate with the plate seal provided. Leave it to stand at 15 ~ 30°C for 2 hrs.
- 9) Washing
Repeat the procedure 3) above to wash the microplate wells.
- 10) Addition of Enzyme substrate
Add 100 μ L each of Enzyme substrate in all wells.
- 11) Enzyme reaction
Cover the microplate with the plate seal provided. Leave the microplate to stand at 15 ~ 30°C in the dark for 30 min.
- 12) Addition of the Reaction stopper
Remove the plate seal and add 100 μ L each of Reaction stopper in all wells and thoroughly mix.
- 13) Absorbance measurement
Measure absorbance of each well by a microplate reader (wavelength 450 nm). When a dual wavelength microplate reader is used, set the reference wavelength at 620 nm or longer. Absorbance must be measured within 30 min after stopping the enzyme reaction.

VI. Determination

1. Preparation of the calibration curve by Mouse HGF standard solution

- 1) Calculate (mean absorbance of the respective concentration of Mouse HGF standard solution - mean absorbance of 0 ng/mL of Mouse HGF standard solution) (Net OD).
- 2) Plot on logarithmic scales Net OD of the respective concentration.

2. Calculation of the mouse HGF in samples

- 1) When absorbance of samples exceeds that of maximum concentration of Mouse HGF standard solution, dilute samples and repeat the assay again.
- 2) Calculate (mean absorbance of samples - mean absorbance of 0 ng/mL concentration of Mouse HGF standard solution) (Net OD).
- 3) Apply the values calculated in 2) above to the calibration curve and determine mouse HGF concentration of samples.

VII. Warnings and precautions

This kit must be used according to the instructions and for the purpose described in this manual. No result is guaranteed in any use or for any purpose other than those described in this manual.

1. General precautions

- 1) Make sure to return the kit to the 15 ~ 30°C before use.
- 2) Do not mix up kit components of different production lots. (Mouse HGF standard solution carries a different lot number from other components in the kit.)
- 3) Assay strictly as instructed.
- 4) Do not use expired reagents.
- 5) Avoid contamination of the kit reagents with microorganisms.
- 6) Thoroughly wash equipment used for the assay and rinse them with purified water.
- 7) Replace micropipette tips for each sample and reagent.

2. Operational precautions

- 1) Measure absorbance of Mouse HGF standard solution for each assay.
- 2) Once assay is started, complete it within the prescribed time and allow the same length of reaction time for Mouse HGF standard solution and samples.
- 3) Make sure that all reactions take place at 15 ~ 30°C.
- 4) Measure absorbance within 30 min after stopping the enzyme reaction.
- 5) Do not scrape the microplate or touch the bottom of wells. Do not dry the surface of the wells during assay.

3. Handling precautions

- 1) Avoid contact of reagents. If they contact skin, wash with plenty of water. Get medical care if need.
- 2) When discarding Sample diluent and Reaction stopper, run sufficient volume of tap water, as Sample diluent contains 0.1% sodium azide and Reaction stopper is acidic.

- 3) Dispose of container and unused contents in accordance with federal, state and local requirements.

VIII. Storage and shelf life

Store the kit at 2 ~ 10°C in the dark and avoid freezing. This kit is stable for 1 year after the date of manufacture. Validity of kit is shown in the package.

IX. Package

1 kit for 96 tests Code No. 1Z85

Assay Procedures and Well Arrangement

Method A - High salt concentration samples

Well Arrangement		Reagents	HGF standard 1A - 1H, 2A - 2B	Samples 2C ~ 12H
1	Addition of the HGF standard solutions and samples	Mouse HGF standard solution Sample* Sample diluent L Sample diluent H	50 μ L -- -- 50 μ L	-- 50 μ L 50 μ L --
2	1st reaction	20 hrs at 15 ~ 30°C		
3	Washing	5 times		
4	Addition of Anti-HGF rabbit antibody	Anti-HGF rabbit antibody	100 μ L	100 μ L
5	2nd reaction	2 hrs at 15 ~ 30°C		
6	Washing	5 times		
7	Addition of Enzyme labeled antibody	Enzyme labeled antibody	100 μ L	100 μ L
8	Third reaction	2 hrs at 15 ~ 30°C		
9	Washing	5 times		
10	Addition of Enzyme substrate	Enzyme substrate	100 μ L	100 μ L
11	Enzyme reaction	30 min in the dark at 15 ~ 30°C		
12	Addition of Reaction stopper	Reaction stopper	100 μ L	100 μ L
13	Absorbance measurement	450 nm		
14	Determination			

*Sample or sample diluted with Sample diluent H.

Method B - Low salt concentration samples

Well Arrangement		Reagents	HGF standard 1A - 1H, 2A - 2B	Samples 2C ~ 12H
1	Addition of the HGF standard solutions and samples	Mouse HGF standard solution Sample* Sample diluent L	50 μ L -- 50 μ L	-- 50 μ L 50 μ L
2	1st reaction	20 hrs at 15 ~ 30°C		
3	Washing	5 times		
4	Addition of Anti-HGF rabbit antibody	Anti-HGF rabbit antibody	100 μ L	100 μ L
5	2nd reaction	2 hrs at 15 ~ 30°C		
6	Washing	5 times		
7	Addition of Enzyme labeled antibody	Enzyme labeled antibody	100 μ L	100 μ L
8	Third reaction	2 hrs at 15 ~ 30°C		
9	Washing	5 times		
10	Addition of Enzyme substrate	Enzyme substrate	100 μ L	100 μ L
11	Enzyme reaction	30 min in the dark at 15 ~ 30°C		
12	Addition of Reaction stopper	Reaction stopper	100 μ L	100 μ L
13	Absorbance measurement	450 nm		
14	Determination			

*Sample or sample diluted with Sample diluent L.

Manufacturer and sold by;

Institute of Immunology Co., Ltd.

1-1-10, Koraku, Bunkyo-Ku,

Tokyo 112-0004, JAPAN

Tel: +81-3-3814-4081

Fax: +81-3-3814-5957