

仕 様 書

ヘリコバクター・ピロリ抗原-JP

コード	6Z32
内容	ウレアーゼ A サブユニット、ウレアーゼ B サブユニット、HSP を含む
由来	菌株 TMPA25 を培養して得られた菌体
精製	超音波破碎および遠心処理
溶媒	Saline (0.1% NaN ₃ 含有)
蛋白濃度	1.0 mg/mL
容量	1.0 mL
保管	-70 °C以下（凍結融解の繰り返しは避けてください）
用途	ELISA、Western blotting

備考 ピロリ菌は超音波により不活化されています。



株式
会社

特殊免疫研究所

〒112-0004 東京都文京区後楽 1-1-10

Tel: 03-3814-4081 Fax: 03-3814-5957

e-mail: info@tokumen.co.jp

<https://www.tokumen.co.jp/>

本製品は研究用にのみ使用できます。人やその他の動物の診断・治療に使用することはできません。

6Z322505-04

Specification Sheet

Helicobacter Pylori Antigen-JP

Code	6Z32
Description	This antigen contains Urease A Subunit, Urease B Subunit, and HSP.
Source	Bacteria obtained by culturing strain TMPA25
Purification	Ultrasonication and centrifugation
Form	Saline (containing 0.1% NaN ₃)
Protein Concentration	1.0 mg/mL
Package	1.0 mL
Storage	Below -70°C. Avoid repeated freeze/thaw cycles.
Application	ELISA, Western blotting

Note The helicobacter pylori is inactivated by ultrasonication.



INSTITUTE OF IMMUNOLOGY Co., LTD.

1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004

Tel: +81-3-3814-4081 Fax: +81-3-3814-5957

e-mail: info@tokumen.co.jp

<https://www.tokumen.co.jp/en>

This product is for in vitro research use only, and is not intended for use in the diagnosis or treatment of humans or animals.

6Z322505-04